

Table S1. Primers and probes used in this study for *Mycoplasma conjunctivae* detection and Sanger DNA sequence analysis.

Primers and probes	Sequences (5'-3')	Use
LPPS-TM-L	CAGCTGGTGTAGCACTTTTTGC	qPCR
LPPS-TM-R	TTAACACCTATGCTCTCGTCTTTGA	qPCR
LPPS-TM-FT	TGCTTCGACTACCAAATATGATGGTGATCCTCT ^a	qPCR - probe
Serstart3	TTTAGTAGACTCCACTTCACC	PCR
Serstart2	CACTATACTTAACAGATAGTCC	Nested PCR, Sequencing
Serstart0	ATACTCAAAGTGGAAATAATGGAA	Sequencing
Serend0	GCAACAACAATAGTAAGAGCAG	Sequencing
lppTA2	TTTGATCTCTCCACCTTCAGC	PCR
lppTA	GGCACTAATAGTGCCTAATTC	Nested PCR

^a5' 6FAM reporter dye and 3' TAMRA quencher.